Intratympanic Needle

Instructions for Use

EN	Intratympanic Needle	Instructions for Use
SV	Intratympanisk Nål	Bruksanvisning
DA	Intratympanisk Nål	Brugsanvisning
FI	Intratympaninen Neula	Käyttöohjeet
NO	Intratympanisk Nål	Bruksanvisning
DE	Intratympanische Nadel	Gebrauchsanweisung
NL	Intratympanische Naald	Gebruiksaanwijzing
IT	Ago Intratimpanico	Istruzioni per l'uso
FR	Aiguille Intratympanique	Instructions d'utilisation
ES	Aguja Intratimpánica	Instrucciones de Uso
PT	Agulha Intratimpânica	Instruções de Uso
PL	Igła Intratympaniczna	Instrukcja użytkowania
CS	Intratympanická jehla	Návod k použití

Medioplast AB
Bronsåldersgatan 2
SE-213 76 Malmö, Sweden.
Phone: +46 40 671 23 00
www.medioplast.com

EN

Intended Purpose

The Intratympanic Needle is a sterile, single use, surgically invasive device intended to penetrate the tympanic membrane and facilitate an intratympanic injection.

Intended Medical Indication

The injection with steroids or antibiotics is typically performed to manage certain inner ear pathologies. Indications for intratympanic injections include sudden sensorineural hearing loss, Ménière's disease, autoimmune inner ear disease and acoustic trauma.

Intended users and required training

Intratympanic administrations are performed by medical doctors, with specialist training in ear, nose and throat disorders or, in some countries, audiological physicians.

Intended Patient Population

Intended patient population is adults with diagnosed inner ear pathology.


Product Description

The tip of the intratympanic needle is sharpened like a standard needle (Quincke type of sharpening) in order to penetrate the tympanic membrane. The intratympanic needle is, after penetrating the tympanic membrane, used to administer a drug to the middle ear. Intratympanic administration is a specific form of drug delivery that is only performed by an ENT specialist via the tympanic membrane.

Technical specifications

The safety stop is designed to indicate when the surface of the intratympanic membrane is reached.

Catalogue No.	Inner Ø	Outer Ø	Length
6066590	0,4 mm	0,7 mm	90 mm

Contraindications

- Uncooperative patient: absolute patient cooperation is required to ensure safe and successful administration.
- A patient with an intact tympanostomy tube (corticosteroid can be instilled via the tube).
- Acute otitis externa.
- An intratympanic tumor or vascular abnormality.
- Suboptimal view of the important landmarks of the tympanic membrane.

Warnings

- DO NOT REUSE and DO NOT RESTERILIZE by any method. The device is a single-use device. If reused or resterilised Medioplast cannot guarantee the functionality nor the safety of the product. Reuse may lead to infection.
- DO NOT push the needle past the safety stop as it may cause permanent damage to vital parts in the middle ear.
- DO NOT use the product after the expiration date. Use after passed expiration date may lead to infection.
- DO NOT use the product if the sterile packaging has been damaged or unintentionally opened before use. Unsterile product may cause infection.

Instructions for use

Preparation

- Before starting the procedure, ensure that the use-by date has not expired and that the sterile package IS NOT damaged
- Prepare and handle the pharmaceutical in accordance with applicable instructions for use.
- Position the patient in a comfortable supine position, with their head rotated 45 degrees to the contralateral ear.
- Use a binocular otology microscope and speculum to visualize the anteroinferior quadrant of the tympanic membrane.
- Wax and debris in the external auditory canal may need to be removed to ensure good visualization of the entire tympanic membrane.
- Apply a suitable anaesthetic.
- Apply a suitable syringe with steroids or antibiotics to the Luer hub.

Procedure Step-by-Step

- Place the Intratympanic Needle, visually confirm that the needle IS NOT inserted further than the safety stop. See Figure 1 Tip of the Intratympanic needle.
- Inject the pharmaceutical using the Intratympanic Needle through the anteroinferior anesthetized region into the middle ear space until full.
- Consider making a secondary perforation, known as a borehole, made on the tympanic membrane in order to relieve pressure, and to avoid barotrauma to the round window.
- Alternatively, a pressure-equalization tube (grommet) can be placed for patients requiring regular intratympanic injections.
- Following administration, the needle is removed, and the patient remains supine with the treated ear facing upwards. Avoid swallowing, yawning, or speaking for

30-40 minutes to facilitate pharmaceutical passage across the round window membrane and prevent leakage into the eustachian tube.

Complications

The most common side effects are transient dizziness, injection site pain, and a burning sensation. Other possible complications or side effects include: Pain, Ear fullness, Vertigo (generally temporary), Headache, Dizziness, Persistent tympanic membrane perforation, Tinnitus, Infection, Syncopal episode, Hearing loss, Tongue numbness

Product disposal

Always follow medical practice and national requirements regarding sharp objects and biohazard when disposing of a used medical device.

Accessories

- Appropriately sized ear speculum
- Binocular otology microscope
- Syringe with standard Luer connection filled with pharmaceuticals.

Reporting

In case of a serious incident in relation to the device, report immediately to Mediplast and the local competent authority.

Note: Duplicate copies of this instruction for use can be ordered free of charge.
The IFU is available for download on the Mediplast website.

This instruction for use has been prepared in English and subsequently translated into other languages. In the event of inconsistencies of interpretation, the English version shall be the definitive version.

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